

CLEANING LABORATORY EVALUATION SUMMARY

SCL #: 2001

DateRun: 06/06/2001

Experimenters: Jason Marshall

ClientType: Electronics Manufacturer

ProjectNumber: Project #1

Substrates: Alloys, Aluminum

PartType: Coupon

Contaminants: Fluxes, Dirt

Cleaning Methods: Ultrasonics

Analytical Methods: Visual, microscopic

Purpose: To evaluate selected cleaner at a lower concentration.

Experimental Procedure: One cleaner, selected by the client based on the previous trials was diluted to 1% using DI water and placed in a VWR Scientific Products Aquasonic 150 HT ultrasonic tank and degassed for 5 minutes. One part was analyzed prior to cleaning to determine a baseline level. The part was then cleaned in this solution for 5 minutes at room temperature. Following cleaning, the part was rinsed for 2 minutes using a Crest 40 kHz ultrasonic tank model 4Ht 1014-6 with DI water at room temperature and then dried for 30 minutes at 160 F. After the part cools to room temperature, a final microscopy analysis was performed.

Results: After analysis, the part looked clean along the seam. There was no major alteration of the surface. Pictures were taken of the part before and after.

Summary:

Substrates:	Alloys, Aluminum				
Contaminants:	Fluxes, Dirt				
Company Name:	Product Name:	Conc.:	Efficiency:	Effective:	Observations:
Tower Products Inc	Tower 270 Cleaner Concentrate	1	0.00	☒	

Conclusion: Cleaning the supplied part with Tower 270 at 1%, 68 F for 5 minutes appears to have been effective.