

# CLEANING LABORATORY EVALUATION SUMMARY

SCL #: 2001  
 DateRun: 11/07/2001  
 Experimenters: Jason Marshall  
 ClientType: Metal  
 ProjectNumber: Project #1  
 Substrates: Steel  
 PartType: Coupon  
 Contaminants: Coatings  
 Cleaning Methods: Ultrasonics  
 Analytical Methods: Gravimetric

Purpose: To compare ultrasonic and immersion cleaning.

Experimental Procedure: Two solutions for each of the five cleaners were diluted to 5% using DI water in 600 mL beakers. One set of beakers were heated to 130 F in a Crest 40 kHz ultrasonic tank. The second set was heated to 130 on a hot plate. Thirty preweighed coupons were coated with the Houghton Veto C3 contaminant and allowed to sit for 2 hours. A second weighing was performed. Three coupons were cleaned in each solution for 3 minutes. Fifteen coupons were cleaned in the ultrasonic tank and the other coupons were cleaned using stir-bar agitated immersion. After cleaning, coupons were rinsed in tap water at 120 F for 30 seconds and dried using a heat gun at 500 F for 1 minute. Once the coupons returned to room temperature, final weights were measured and efficiencies were calculated.

Results: The ultrasonic energy was found to be more successful than the immersion cleaning was. The following table lists the calculated results for both methods.

Cleaner	Contaminant	Coupon 1	Coupon 2	Coupon 3	Average	Previous Results
Multikleen	Ultrasonics	99.99	99.77	99.83	99.86	97.88
Inproclean	Ultrasonics	99.93	99.56	100.01	99.83	99.32
Beyond	Ultrasonics	99.87	100.11	100.05	100.01	99.26
Dasco	Ultrasonics	100.02	100.03	100.04	100.03	93.11
Certa	Ultrasonics	100.06	99.70	100.03	99.93	96.86
Multikleen	Immersion	99.88	99.65	99.17	99.57	99.59
Inproclean	Immersion	99.93	99.70	99.58	99.74	99.99
Beyond	Immersion	99.71	98.72	99.48	99.30	99.44
Dasco	Immersion	98.74	97.53	98.22	98.16	99.27
Certa	Immersion	99.09	99.26	99.64	99.33	98.15

Comparison of this data to the previous trial data revealed that the ultrasonic cleaning was more effective in a cleaner solution. The new immersion cleaning results were nearly identical to the previous trial results.

Summary:

<b>Substrates:</b>	Steel				
<b>Contaminants:</b>	Coatings				
Company Name:	Product Name:	Conc.:	Efficiency:	Effective:	Observations:
Heatbath Corporation	Multi-Kleen 1568	5	99.86	<input checked="" type="checkbox"/>	ultrasonics
Heatbath Corporation	Multi-Kleen 1568	5	99.57	<input type="checkbox"/>	Immersion
Oakite Products	Inproclean 3800	5	99.83	<input checked="" type="checkbox"/>	ultrasonics
Oakite Products	Inproclean 3800	5	99.74	<input type="checkbox"/>	Immersion
Today & Beyond	Beyond 2001	5	100.01	<input checked="" type="checkbox"/>	ultrasonics
Today & Beyond	Beyond 2001	5	99.30	<input type="checkbox"/>	Immersion
DA Stuart Company	Dasco Kleen 3250	5	100.03	<input checked="" type="checkbox"/>	ultrasonics
DA Stuart Company	Dasco Kleen 3250	5	98.16	<input type="checkbox"/>	Immersion
Houghton International	Cerfa Kleen 5387	5	99.93	<input checked="" type="checkbox"/>	ultrasonics
Houghton International	Cerfa Kleen 5387	5	99.33	<input type="checkbox"/>	immersion

Conclusion:

## **CLEANING LABORATORY EVALUATION SUMMARY**

When comparing similar bath solutions (fresh), ultrasonic cleaning was found to be more effective than the immersion cleaning was for removing the rust preventative.