

CLEANING LABORATORY EVALUATION SUMMARY

SCL #: 2002
 DateRun: 04/29/2002
 Experimenters: Jason Marshall
 ClientType: Metal
 ProjectNumber: Project #1
 Substrates: Aluminum
 PartType: Coupon
 Contaminants: Adhesive
 Cleaning Methods: Immersion/Soak
 Analytical Methods: Gravimetric
 Purpose: To identify possible alternatives to acetone for adhesive removal

Experimental Procedure: Nine products were selected from the lab's database of alternative cleaning products based on client input. Six semi-aqueous and three aqueous products were chosen for evaluation. An additional product was supplied by the client. The aqueous products were diluted to 10% using DI water in 600 ml beakers. All ten products were heated to 130 F on a hot plate. Thirty preweighed coupons were coated with Bob Smith Industries Insta-Cure adhesive (7085-85-0) and then sprayed with the Insta-Set (64742-89-8, 99-97-8) and allowed to sit for 10 minutes at room temperature. After drying, the coupons were weighed a second time. Three coupons were placed in each beaker and allowed to sit for 10 minutes, then rinsed in tap water at 120 F for 10 seconds and dried with a Master Appliance heat gun at 500 F for 10 seconds. A third weighing was taken and efficiencies were calculated. The same coupons were placed back into the cleaning products and soaked for an additional 50 minutes. At the end of the hour, coupons were rinsed, dried and weighed again. Final efficiencies were recorded.

Results: The client supplied product and DuPont's DBE performed similarly, causing the adhesive to become soft and liquefied. Both products are ester based. The International Products Surface Cleanse 930 was found to cause the adhesive to lift off the metal surface in one big stiff piece. The US Polychemical Co Polyspray Jet 790 XS had similar results after the 60 minutes of cleaning. Two products, Westech 1050 and RST-5 attacked the aluminum coupons.

The 10 minute cleaning results are listed in Table 1.

Cleaner	Coupon 1	Coupon 2	Coupon 3	Average	Std Dev	Observations
Bio T Max	-0.33	3.79	1.65	1.70	2.06	
Westec 1050	-93.88	-77.00	-75.24	-82.04	10.29	Attacked aluminum
DBE	-109.73	-174.95	-153.35	-146.01	33.22	sticky/soft
RST-5	-1.48	20.18	15.56	11.42	11.41	
Ionox HC2	28.53	18.93	24.39	23.95	4.81	
Shopmaster LpH	-1.83	1.06	4.30	1.18	3.07	
Surface Cleanse	0.47	-10.35	-17.15	-9.01	8.89	wipe/peel off
Polyspray	41.97	24.67	2.84	23.16	19.61	wipe/peel off
Santosol DME	-146.75	-130.63	-202.49	-159.96	37.71	sticky/soft
Bio T Max	0.61	3.58	0.67	1.62	1.69	

After cleaning for an hour, little improvement was observed for most of the products.

Table 2. Cleaning Efficiencies @ 60 Minutes

Cleaner	Coupon 1	Coupon 2	Coupon 3	Average	Std Dev	Observation
Bio T 200 A	0.04	3.35	1.51	1.63	1.66	
DBE	-237.78	-324.05	-266.71	-276.18	43.90	
RST-5	21.37	109.68	114.16	81.74	52.33	2/3 peeled but attacked metal
Ionox HC2	43.03	37.41	48.90	43.11	5.75	

CLEANING LABORATORY EVALUATION SUMMARY

Shopmaster LpH	0.56	1.87	3.55	1.99	1.50	
Surface Cleanse	92.67	97.18	94.40	94.75	2.27	3/3 peeled
Polyspray	46.10	87.50	95.50	76.37	26.52	2/3 peeled
Santosol DME	-216.06	-203.42	-338.42	-252.63	74.56	
Bio T Max	0.25	4.48	1.55	2.09	2.17	

Summary:

Substrates:		Aluminum			
Contaminants:		Adhesive			
Company Name:	Product Name:	Conc.:	Efficiency:	Effective:	Observations:
Bio Chem Systems	Bio T 200 A	100	1.63	<input type="checkbox"/>	
Chem Tech Solutions	1050 Magic	100	-82.04	<input type="checkbox"/>	attacked aluminum
Invista S.a.r.l	Flexisolv DBE Ester	100	-276.18	<input type="checkbox"/>	sticky/soft
Universal RenigIngsmitel - Mulder Hardenberg	RST 5	100	81.74	<input type="checkbox"/>	Attacked metal
Kyzen Corporation	Ionox HC 2	100	43.11	<input type="checkbox"/>	
Buckeye International	Shopmaster LPH	10	1.99	<input type="checkbox"/>	
International Products Corporation	Surface Cleanse Concentrated Neutral 930	10	94.75	<input checked="" type="checkbox"/>	adhesive peeled off
US Polychem Corporation	Polyspray Jet 790 XS	10	76.37	<input type="checkbox"/>	Two coupons had adhesive peel off after 60 minutes.
Bio Chem Systems	Bio T Max	100	2.09	<input type="checkbox"/>	
Surface Specialties Inc	Santosol DME	100	-252.63	<input type="checkbox"/>	soft/sticky

Conclusion:

Surface Cleanse 930 will be used in attempt to remove the film from the metal part.